



Direct Mobility Solutions

Manual Wheelchairs



Manual Wheelchairs

Catalyst 4

The Catalyst 4 & 4C are the lightest most durable and energy efficient folders on the market today.

- The Catalyst 4 series axle is mounted less than 1/2" forward to traditional axle positions. The simple axle plate design offers vertical adjustments of the axle in five positions.
- Catalyst 4 & 4C offer the option of a swing away front end featuring our patented footrests or 70° or 80° fixed frames. The swing away hangers come in 4 bends (60°, 70°, 80° & 90°) are light and easy to use and so strong we offer a lifetime warranty.
- Ki standard spoke uses 30 radial spokes to provide a stiff and durable wheel.
- Exclusive to Ki, the magnesium Maxx Mag wheel creates a sleek, modernized look while providing greater stiffness and lower weight than classic nylon composite mag wheels. Maxx Mag wheels are available in 20", 22" and 24".

Starting at \$1,700



Ki Mobility
by Etac

Tracer EX2

Dual-axle allows for repositioning for hemi seat-to-floor heights. Combining a heavy-duty seat inner liner, low-maintenance carbon steel frame and interchangeable components with the Invacare 9000 Wheelchair series, this wheelchair sets the benchmark for long-lasting, durable wheelchairs.

- Attractive, silver vein finish
- Heavy-duty inner liners keep seat and back from stretching
- 14-gauge crossbraces add durability and strength
- Folds to 13" wide for convenient travel and storage
- Dual-embossed upholstery for added reinforcement



Starting at \$658.57



Tracer IV Heavy-Duty

Combining a reinforced, heavy-duty carbon steel frame which includes gusseted sides, caster journals and dual crossbraces, the Invacare Tracer IV Wheelchair is long-lasting and low-maintenance. Heavy-duty seat liner provides added durability. The T4 is customizable as a builder wheelchair to meet the needs of individual clients.

- Heavy-duty inner liners keep seat and back from stretching
- Attractive, silver vein finish
- Padded arms for added comfort
- 14-gauge crossbraces add durability and strength
- Folds to 13" wide for convenient travel and storage



Starting at \$1,568



Manual Wheelchairs

Strongback 24HD (Heavy Duty)

The lightweight and ergonomic design provides exceptional support and improves posture. The STRONGBACK 24HD Heavy Duty is the first lightweight wheelchair with an impressive 375 lb. weight capacity.

- Lightweight frame with an impressive 375lb. (170kg) weight capacity
- Award winning, patented backrest for unmatched postural support
- Ultra comfortable, easy to propel design
- Unmatched attendant control with integrated drum brakes for superior maneuverability

Starting at **\$1,199**



Strongback 24



Experience comfort and mobility with the STRONGBACK Comfort : 24 Wheelchair from Strongback Mobility. Its lightweight and ergonomic design provides exceptional support and improves posture.

- Patented STRONGBACK ergonomics with easily adjustable posture support in the backrest
- Quick-release rear wheels
- Compact foldable design
- Ultra-lightweight
- Swing-away detachable footrests
- Weight capacity: 300 lbs.

Starting at **\$909**



Feather Chair

Feather

At only 13.5 pounds*, the Feather Chair is one of the lightest wheelchairs on the planet. Taking advantage of brand-new materials, the Feather Chair sets a new standard. Practically everyone can lift this wheelchair into their car, truck or SUV. It is super easy to propel and push.

Standard features include wheel locks - easily accessible by the seated user - at the front of the frame. Wheel locks prevent the wheelchair from rolling forward or backward.

For the caregiver, there is an integrated braking system, operated by bicycle style locking brake levers located at the push handles. These handbrakes provide the caregiver control, security, and convenience when descending or pausing on a slope.

Starting at **\$877**



Transport Wheelchairs

Excursion 8



Experience exceptional comfort and convenience with the STRONGBACK Excursion : 8 Transport Chair from Strongback Mobility. Its lightweight design and ergonomic features make it the perfect choice for enhanced mobility.

- Easily adjustable posture support in the backrest
- Compact foldable design
- Ultra-lightweight
- Swing-away detachable footrests
- Ergonomic comfort-grip handles
- Comfortable desktop-length armrests
- Weight capacity: 300 lbs.

Starting at \$779



Feather Transport



At only 13 pounds*, the Feather Transport is the lightest transport wheelchair ever created.

Transporting individuals has never been easier! With brand-new materials, the Feather Transport sets a new standard for transport wheelchairs.

Anyone can lift this wheelchair into a car, truck or SUV.

It is super easy to propel and push.

The chair also folds down into a compact package.

Starting at \$499



Lightweight Transport Chair

with Hand Brakes - 19" with Swing Away Footrests - Blue



The Lightweight Transport Chair with Hand Brakes has a lightweight aluminum frame and padded upholstery. It is foldable and has a quick release fold-down back, as well as padded, fixed full-length arms for added comfort.

- Ideal for transport, travel & safety
- Most versatile design for all terrain
- Comfortable padded arms
- Easy to use footrests
- Weight capacity: 300 lbs.

Starting at \$334.95



Lightweight Transport Chair

with Locking Hand Brakes, Flip Up Arms and Anti-Tippers - Red



We love the 350 Transport Chair because it's lightweight, yet versatile and all-terrain, with locking handbrakes, flip up desk arms, 8" wheels and anti-tippers.

- Lightweight aluminum frame
- Flip-up desk arms
- Hand brakes & locking rear wheels
- Non-ski foot plates and heel loops
- Foldable
- 8" wheels
- Anti-tippers

Starting at \$334.95



Ultra Lightweight Folding Wheelchairs

Catalyst 5



Ki Mobility's flagship product, the Catalyst 5 has since redefined the ultra-lightweight folding market by raising the bar in weight, performance, durability and functionality. Weighing in at 23 lbs and constructed from ultra-lightweight 7000 series aluminum, the Catalyst 5 is arguably the best choice in high-performance folding wheelchairs.

Starting at **\$2,352**

TiLite Aero X



The Aero X features TiLite's Tru-Fit Component System making components clamp directly onto the side frames which eliminating one of the main stress points by reducing the number of holes in the frame by 50%. In addition, the Tru-Fit System allows for infinite adjustability of its components. This design innovation also greatly improves the durability of the frame.

Starting at **\$2,586**

Helio C2



This flagship model is one of the lightest, most technologically advanced wheelchairs in its category. The HELIO C2 is the result of a series of innovations, it features revolutionary components, an impeccable finish and all the advantages carbon fiber offers.

Starting at **\$4,535**

Veloce



Offering unrivaled agility, lightness, stability and responsiveness, this folding chair promises both performance and comfort. Featuring the all-important benefits of a compact folding system, VELOCE is easy to transport and to adjust. In short, it adapts to your active lifestyle.

Starting at **\$5,795**

Ultra Lightweight Rigid Wheelchairs

Rogue 2



Say goodbye to the ordinary and embrace the extraordinary as you explore limitless color combinations to create a wheelchair that truly reflects your style and personality. Whether you're drawn to bold, vibrant hues or prefer sleek, subtle tones, the Color Visualizer brings your vision to life in real time.

Starting at **\$2,873**

Apex C



The APEX Carbon wheelchair by Motion Composites belongs to a new era of rigid wheelchairs. A balanced blend of technology and design, it is the ultimate in lightness, durability, flexibility and style. An high-end technology and all the advantages of carbon fiber in a modern-looking, fully adjustable chair.

Starting at **\$2,595**

TiLite Aero T



Our most advanced aluminum chair, the Aero T's dual-tube frame minimizes weight while maximizing performance. Add to that full adjustability and TiFit, and you've got something special.

Starting at **\$3,191**

TiLite Aero Z



The open frame of the aluminum Aero Z calls out to the open road. The Aero Z is the only TiFit, mono-tube aluminum frame on the market. In a class by itself, the Aero Z is lightweight and fully adjustable and it won't break the bank!

Starting at **\$2,873**



Direct Mobility Solutions

Direct Mobility Solutions is locally owned and operated in the Philadelphia area. We pride ourselves on quality equipment and unmatched service and expertise.

Call Us Today and we will help you choose the best option to fit your lifestyle.
