



Direct Mobility Solutions

Patient Lifts



PL400H



With a maximum weight capacity of 400 lb., the PL400H is easy to handle and maintain while providing robust healthcare options. The hydraulic pump provides a low-cost lifting option while allowing caregivers to easily transfer patients.

- 400 Lb capacity
- Six-Point spreader bar
- Manual Open/Close base legs
- Optional Upgrade to battery power

Starting at \$626

PL400HE



Featuring a 400 lb. maximum weight capacity and Bestcare battery powered control system for power lifting, the PL400HE is easily adaptable for a wide range of challenges while still maintaining a high standard of comfort for patients.

- 400 lb. carrying capacity
- Six-point spreader bar
- Electronic control system
- Optional digital weigh scale
- Optional two-point spreader bar

Starting at \$1,250

PL182



The BestLift® PL182 – Full Body Patient Lift

A perfect blend of form and function for patient transfers. With an overall base length of 44.9", maximum width of 42.9", and base height of 4.5", it is the perfect solution for transferring bariatric patients in small rooms with low beds.

- 400 Lb capacity
- Dual front castors
- New electronics with 24/7 capability
- New hand operated base opening

Starting at \$1,297

PL228 (HD Option)



The BestLift® PL228 – Full Body Patient Lift

A perfect blend of form and function for patient transfers. With an overall base length of 44.9", maximum width of 42.9", and base height of 4.5", it is the perfect solution for transferring bariatric patients in small rooms with low beds.

- 500 Lb capacity
- Dual front castors
- New electronics with 24/7 capability
- New hand operated base opening

Starting at \$1,439



PL400

One of the latest ProCare products. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.

- 400lb capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Optional Scale

Starting at \$1,529



PL500 (HD Option)

One of the latest ProCare products. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.

- 500lb capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Optional Scale

Starting at \$1,982



PL600 (HD Option)

The heaviest duty ProCare patient lift. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.

- 600lb capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Optional Scale

Starting at \$2,709

STA400



The STA400 is Bestcare's entry level solution for Standing Transfer Aids.

It carries a maximum weight of 400 lb. while allowing the user to actively participate in the transfer. Swivel seats supports the patient and allows for a quick and comfortable change from sitting to standing and vice versa. Additional attachment hooks for a support strap are also included.

- 400 lb. carrying capacity
- Folding seats
- Dual Knee pads
- Non-slip foot rest
- Grab bar for user
- Optional Raising Belt for stand transfer and support
- Optional support strap
- Optional Seat lock

Starting at \$689

STA450



The STA450 is Bestcare's premier solution for Standing Transfer Aids.

In addition to the benefits of the STA400, the STA450 has a higher weight capacity at 450 lb. and has an expandable base. Together these two features allow caregivers to provide care for a wider range of residents as well as greater access to residents in more challenging environments.

- 450 lb. carrying capacity
- Manual base open/close
- Folding seats
- Dual Knee pads
- Non-slip foot rest
- Grab bar for user
- Hooks for optional side support strap
- Optional "S" short model featuring a low seat

Starting at \$882



PL400EF

The BestLift™ PL400EF represents a new standard in full body lifts for individuals who travel or have storage issues in their home. It is a unique lifter with a small footprint and foldable design that assembles or folds in 3 simple steps allowing for easy storage or transport. As with all Bestcare lifts, this new Bestcare lifter continues to insure maximum safety and comfort.

- Foldable / Portable / Transportable
- 400lb capacity
- Smaller and lighter
- 6-point spreader bar (2 point optional)
- Fits to car trunk
- New Bestcare Control system

Starting at \$1,342



Sit to Stand Lifts (Hydraulic and Electric)



SA182H/SA228H



The BestStand SA182/SA182H (400lb series) and SA228/SA228H (500lb series) are affordable Care sit-to-stand lift and transfer solutions. They are compact, easy to maneuver, and transportable.

- SA182/SA182H are rated at 400lb safe working load
- SA228/SA228H are rated at 500lb safe working load
- SA182H and SA228H are hydraulic models with optional electric upgrade
- SA182 and SA228 are electric models

Starting at \$1,184

SA182/SA228



The BestStand SA182/SA182H (400lb series) and SA228/SA228H (500lb series) are affordable Care sit-to-stand lift and transfer solutions. They are compact, easy to maneuver, and transportable.

- SA182/SA182H are rated at 400lb safe working load
- SA228/SA228H are rated at 500lb safe working load
- SA182H and SA228H are hydraulic models with optional electric upgrade
- SA182 and SA228 are electric models

(SA182) Starting at \$1,797

(SA228) Starting at \$2,009

SA400H/SA400HE



The BestStand™ SA400H/HE Mini are our latest innovative design for quick and easy sit-to-stand transfer. They are small and compact. Easy to maneuver in tight space. Available in both hydraulic and electric.

- 400lb Capacity
- Unique T-bar design for quick lifting action
- Hydraulic with optional electric upgrade
- Adjustable knee pad
- Detachable foot plate


(SA400H) Starting at \$1,210

(SA400HE) Starting at \$1,814

Best Patient Sling for You...

Bestcare
L.L.C.

Universal Sling




Universal slings or "U" shaped slings mold around the body easily. Features six to eight hook-up straps for an array of attachment options to a six-point spreader bar. More attachment straps allow adjustments for better weight distribution and more comfortable positions. Padded legs reduce skin injury and discomfort during sling application. Select Mesh models for wet applications and quicker drying times. Mesh models come without padded legs. Universal slings are available with or without head support.

- Head Support has 8 Straps / No Head Support has 6 Straps
- Compatible with 6Pt. / 4Pt. / 2Pt. Spreader Bar
- Fits on Full Body Patient Lift

Fabric: **Polyester** **Mesh**

Pediatric Sling




Pediatric slings are specialty slings designed for use with children. This is a "hammock" style sling with added upper thigh support. Unlike the common split-leg sling, each thigh is supported individually and both thighs are kept close together by an additional secure loop. Other features include padded legs, removable head-shell, security belt and pleasant nursery graphics. The Pediatric sling is also available in Mesh for wet applications.

- Six Straps and one loop
- Compatible with 4Pt. Spreader Bar
- Fits on Full Body Patient Lift

Fabric: **Polyester** **Mesh**

Sani Sling




Sani-Slings are made for maximum functionality and patient dignity. Eliminates need for poorly designed slings with commode opening. Easy to apply and use, these slings are well padded under the arms and lumbar region for maximum comfort. The Sani-Sling also works extremely well on stand assist lifts for patients requiring lower back and body support during a sitting transfer.

- Four Straps
- Compatible with 6Pt. / 4Pt. / 2Pt. Spreader Bar and 4Pt. Boom
- Fits on Full Body Patient Lift and Stand Assist Lift

Fabric: **Polyester**

Replacement Sling for Liko® / Guldmann®




Replace your Liko® or Guldmann® split-leg sling with this compatible sling. These slings are made from high-weight polyester fabric and are available with or without head support.

- Four Straps
- Compatible with 6Pt. / 4Pt. / 2Pt. Spreader Bar
- Fits on Full Body Patient Lift

Fabric: **Polyester** with PU Coating

Replacement Sling for Hoyer®




Replace your Hoyer® classic sling with this compatible sling. These slings feature padded body and legs and are available with or without head support. Head support versions feature removable plastic backed head-shell.

- Four Straps
- Compatible with 6Pt. / 4Pt. / 2Pt. Spreader Bar
- Fits on Full Body Patient Lift

Fabric: **Polyester**

Replacement Sling for Invacare®




Replace your full body Invacare® sling with this compatible sling. These slings are available in Mesh fabric with or without a commode opening. Solid fabric and bariatric sizes are also available.

- Four Straps
- Compatible with 6Pt. / 4Pt. / 2Pt. Spreader Bar
- Fits on Full Body Patient Lift

Fabric: **Polyester** **Mesh**

Stand Assist Sling




Experience the most comfortable stand assist sling on the market! Lumbar and underarm padding provide maximum comfort for the patient. For those that require additional support an optional buttock strap can be used in conjunction with this sling. A calf strap is available for lower leg stabilization.

- Two Straps
- Compatible with 4Pt. Boom
- Fits on Stand Assist Lift

Fabric: **Polyester**

Disposable Sling (Single Patient Specific Sling)



Bestcare Single Patient Slings are made from a polypropylene nonwoven layered fabric which allows for a cost effective method of controlling cross contamination and serves as an excellent alternative to woven cloth slings whose cost accumulates over time along with required upkeep and laundering. When more slings are available, staff can assign a sling to each individual Patient. No more sharing slings. No chasing down slings that have gone missing with the laundry. Latex free and tested at 1.5 times the safe working load. SPS slings are intended to be used with a single patient only and cannot be laundered. When soiled, dispose of the sling accordingly.

- Liko® / Guldmann® Replacement Disposable Sling has 4 Straps
- Invacare® Replacement Disposable Sling has 4 Straps.
- Universal Disposable Sling Head Support has 8 Straps / No Head Support has 6 Straps
- Repositioning Disposable Sling has 10 Straps
- Stand Assist Disposable Sling has 2 Straps
- Compatible with 6Pt. / 4Pt. / 2Pt. Spreader Bar. **Only Stand Assist Disposable Sling compatible with 4Pt. Boom**

Fabric: **Nonwoven Disposable** **Special Woven Disposable**

Universal Sling				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-UM882	✗	Mesh	M	400LB
SL-UM883	✗		L	450LB
SL-UP892	✗	Polyester	M	400LB
SL-UP893	✗		L	450LB
SL-UP894	✗		XXL	600LB
SL-UMH831	✓	Mesh	S	400LB
SL-UMH832	✓		M	400LB
SL-UMH833	✓		L	450LB
SL-UPH861	✓	Polyester	S	400LB
SL-UPH862	✓		M	400LB
SL-UPH863	✓		L	450LB
SL-UPH864	✓		XXL	600LB

Pediatric Sling				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-PM111	✓	Mesh	Standard	100LB or Less
SL-PP121	✓	Polyester	Standard	100LB or Less

Universal Sling				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-TP331	✗	Polyester	S	400LB
SL-TP332	✗		M	400LB
SL-TP333	✗		L	450LB
SL-TP334	✗		XXL	600LB

Stand Assist Sling				
Item No.	Description	Material Type	Size	Weight Capacity
SL-SA661	Deluxe Padded	Polyester	S	400LB
SL-SA662			M	400LB
SL-SA663			L	450LB
SL-SA664			XXL	600LB
SL-SA669	Buttock Support Strap for Stand Assist		Standard	400LB
SL-SA669B			XXL	600LB
SAL1230-KB	Stand Assist Knee Strap			Standard

Disposable Replacement Sling for Liko® / Guldmann®				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-GD431	✗	Nonwoven Disposable	S	600LB
SL-GD432	✗		M	
SL-GD433	✗		L	
SL-GD434	✗		XL	
SL-GDH441	✓	Nonwoven Disposable	S	
SL-GDH442	✓		M	
SL-GDH443	✓		L	
SL-GDH444	✓		XL	


Replacement Sling for Liko® / Guldmann®				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-GP471	✗	Polyester with PU Coating	S	600LB
SL-GP472	✗		M	
SL-GP473	✗		L	
SL-GP474	✗		XXL	
SL-GPH481	✓		S	
SL-GPH482	✓		M	
SL-GPH483	✓		L	
SL-GPH484	✓		XXL	

Replacement Sling for Hoyer®				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-HC70002	✗	Polyester	M	500LB
SL-HC70001	✗		L	500LB
SL-HC70000	✗		XL	600LB
SL-HC70012	✓		M	500LB
SL-HC70011	✓		L	500LB
SL-HC70010	✓		XL	600LB

Replacement Sling for Invacare®					
Item No.	Commode Opening	Material Type	Color Code	Size	Weight Capacity
SL-R112	✗	Polyester	Yellow	M	450LB
SL-R113	✗		Green	L	
SL-R117	✗		Blue	XL	
SL-R110	✗	Mesh	Yellow	M	600LB
SL-R111	✗		Green	L	
SL-R140	✗		Black	XXL	
SL-R114	✓	Mesh	Yellow	M	450LB
SL-R115	✓		Green	L	
SL-R116	✓		Blue	XL	
SL-R141	✓		Black	XXL	
					600LB

Universal Disposable Sling				
Item No.	Head Support	Material Type	Size	Weight Capacity
SL-UD812	✗	Nonwoven Disposable	M	600LB
SL-UD813	✗		L	
SL-UDH821	✓		S	
SL-UDH822	✓		M	
SL-UDH823	✓		L	
SL-UDH824	✓		XXL	

	Small	Medium	Large	X-Large	XX-Large
Color Code	Red	Yellow	Green	Blue	White



Size & Body Weight

	Small	Medium	Large	Bariatric
A-Circumference around the shoulders (Inches)				
Universal Split Leg	24-45	37-57	47-67	57-80
Pediatric	18-28	N/A	N/A	N/A
B-Circumference around the chest (Inches)				
Sani(Toileting)	18-40	16-46	28-50	34-73
Stand Assist	26-42	30-46	34-53	46-64
Body Weight (LB)				
	Small	Medium	Large	Bariatric
Universal Split Leg	85-140	140-200	200-450	400-600
	Medium	Large	X-Large	XX-Large
Replacement of Invacare Sling	100-195	150-300	210-400	290-600



Polyester



NONWOVEN DISPOSABLE

-infection control
-polypropylene
-latex free



Mesh

FOR BATHING APPLICATIONS
-woven fabric that has a large number of closely spaced holes.



SPECIAL WOVEN DISPOSABLE

-polyester
-acrylic coated fabric



Polyester with PU Coating

Disposable Replacement Sling for Invacare®				
Item No.	Commode Opening	Material Type	Size	Weight Capacity
SL-RD114	✓	Special Woven Disposable	M	600LB
SL-RD115	✓		L	
SL-RD116	✓		XL	
SL-RD141	✓		XXL	

Stand Assist Disposable Sling				
Item No.	Description	Material Type	Size	Weight Capacity
SL-DSA662	Deluxe Padded (Toileting)	Nonwoven Disposable	M	400LB
SL-DSA663			L	

Repositioning Disposable Sling			
Item No.	Material Type	Size	Weight Capacity
SL-RPD1000	Special Woven Disposable	Standard	1000LB

1000LB safe work load ONLY when all loop straps are properly attached.



Direct Mobility Solutions

Direct Mobility Solutions is locally owned and operated in the Philadelphia area. We pride ourselves on quality equipment and unmatched service and expertise.

Call Us Today and we will help you choose the best option to fit your lifestyle.
